

3.5	Ini to 2nd	2.0	2.5	2.8	3.3	3.8	4.3		3405	3805	4305	4805
D1	●	●	●	●	◐	-	-		-	●	-	
D2	●	●	●	-	-	-	-		●	●!!	-	
D3	●	●	◐	-	-	-	-		●	-	-	
D4	●	●	◐!!	-	-	-	-		-	-	-	

4.0	Ini to 2nd	2.0	2.5	2.8	3.3	3.8	4.3		3405	3805	4305	4805
D1	●	●	●	●	●	◐	-		-	-	●	
D2	●	●	●	-	◐	-	-		-	●	◐!!	
D3	●	●	●	-	◐	-	-		-	●!!	-	
D4	●	●	●	◐	-	-	-		-	-	-	

4.5	Ini to 2nd	2.0	2.5	2.8	3.3	3.8	4.3		3405	3805	4305	4805
D1	●	●	●	●	●	●	◐		-	-	-	
D2	●	●	●	●	●	◐	-		-	-	●!!	
D3	●	●	●	●	-	◐	-		-	-	●!!	
D4	●	●	●	●	◐!!	-	-		-	-	-	

5.0	Ini to 2nd	2.0	2.5	2.8	3.3	3.8	4.3	4.8	3405	3805	4305	4805	5305
D1	●	●	●	●	●	●		◐					
D2	●	●	●	●	●	●	◐					●!!	
D3	●	●	●	●	●	-	◐					●!!	
D4	●	●	●	●	●	◐!!							

Limitation of Liability

This instruction is a guideline only. The surgeon, therefore, takes full medical responsibility for the application of the surgical drilling protocol.